

COVID19 - RBD Rapid Test



Model/s:

SARS-CoV-2 RBD Protein IgG Detection Kit (Colloidal Gold-Based)

Intended Use:

This product is intended for the in vitro qualitative detection of RBD (SARS-CoV-2 Spike Receptor-Binding Domain) IgG antibody in human clinical samples (Serum/Plasma/Whole Blood). For rapid detection within 15 min. No testing equipment required.

Size: : 20 test / box

Ref. Code : E2602C

Manufacturing : Vazyme Medical Technology Co. Ltd.
Nanjing, China

Storage : 4°C ~ 30°C

Shelf Life : 18 Months

Certifications : CE



PT Fokus Keluarga Sehat

www.fkssehat.com

info@fkssehat.com

+62 813 8005 8788

FKSsehat