



# DENGUE IgG/IgM & NS1 RAPID TEST

(Whole Blood/Serum/Plasma)

- SIMPLE TO USE
- HIGH ACCURACY
- FAST READOUT
- DETECT DRUG ABUSE UP TO 3-5 DAYS

## About Dengue

Dengue viruses are spread to people through the bite of an infected Aedes species mosquito. Dengue is common in more than 100 countries around the world. 40% of the world's population, about 3 billion people, live in areas with a risk of dengue.

## Why do we need Dengue Test?

An accurate diagnosis of Dengue at an early stage is important not only for an Etiological diagnosis for the patient but also to identify individuals that might serve as a source of infection.

## Product Information:



The AyoTes Dengue IgG/IgM & NS1 Rapid Test (Whole Blood/Serum/Plasma) is a rapid chromatographic immunoassay for the qualitative detection of Dengue virus (DENV) NS1 antigen, anti-DENV IgG and anti-Dengue IgM in human whole blood, serum, or plasma as an aid in the diagnosis of primary and secondary Dengue infections.

The designed cutoff concentrations and direct calibration for these drugs are as follows

Cat no.	Specimen	Shelf Life	Storage
IDAT-AS2023005	Serum/Plasma/Whole Blood	24 Months	2-30°C




## How to Use:



## Results interpretation: For Dengue IgG/IgM Test



## Results interpretation: For Dengue NS1



### Important:

Do not interpret the results after 5 minutes

For orders & inquires

**PT Fokus Keluarga Sehat**

Jl. KH. Mansyur No. 129-130, RT.9/RW.3. Karet Tengsin, Tanah Abang District, Central Jakarta City  
Special Capital Region of Jakarta 10250, Indonesia

0813-8005-8788

info@fkssehat.com